

Tuttnauer Clinic

LEADING IN
STERILIZATION
SINCE

1925

T-Top Autoclave

Efficiency with maximum reliability

Featuring the only IoT software designed for the smart clinic, connecting sterilization devices and communicating on-site or remotely.



Designed to Meet Your Healthcare
Clinics' Everyday Sterilization Needs

Affordable and reliable for healthcare clinic sterilization

Tabletop Medical Class B Autoclave From The World's Premium Manufacturer



Efficient



Reliable Performance



Plug & Play



Unique User Interface

T-Top is a family of tabletop autoclaves designed specifically for medical and dental clinics in small- to medium-size practices. Developed with simplicity, efficiency, and cost-effectiveness in mind, **T-Top** devices address the core sterilization needs of everyday clinical workflows, helping increase productivity and streamline routine operations.

Available in multiple chamber configurations, **T-Top** delivers the flexibility and capacity clinics need in an affordable, reliable solution they can depend on for daily sterilization workloads.

Combining performance, ergonomics, advanced technology, and ease of use, **T-Top** offers a well-balanced and robust autoclave platform that meets the requirements of standard medical and dental practices.



Sterilizer with high performance

Complete solution for clinic instrument reprocessing.

Making Your Process Easier, Safer and Affordable.

User friendly advanced technology with Wi-Fi connectivity and intuitive colour touch screens.



Plug and Play

No dedicated installation connections needed, simply a power source

- ◆ Multicolour user-friendly touch screen with intuitive icons
- ◆ Easily start and run programs
- ◆ Compact and ergonomic, designed to fit into your space



Reliable Performance

Quality product that will do all you need and expect

- ◆ Easy to operate, maintain and use
- ◆ Internal memory - storing up to 500 cycles
- ◆ Advanced technology to deliver longevity and reliability



Safety & Regulation

Meets strictest international standards.



Traceability & Connectivity

Advanced Technology that provides real-time data intelligence

- ◆ 2 USB Ports
- ◆ Wi-Fi connection & T-Connect software
- ◆ Monitor devices remotely and download device data
- ◆ Remotely access cycle history, get notifications and alerts.
- ◆ Technicians and service centers can remotely analyze autoclaves, diagnose and resolve problems
- ◆ Paperless documentation
- ◆ R.P.C.R



Innovation

Best in class technology

- ◆ Programming and customization functions
- ◆ Specially designed user interface to speed up and optimize daily workload
- ◆ Time saver setting to set tests and sterilization cycles to start at any time

Designed for your ease of use

A blend of perfect performance and functionality to **make your user experience intuitive and hassle free.**

Sterilization at the Click of a Button

Just select the required program and press "start", a click and go feature to help speed up and optimize your instrument reprocessing workload.



T-Top simplifying your sterilization process



Documentation and Traceability



- ◆ Advanced documentation and traceability options, including barcode
- ◆ Personalized username and password



Printouts for saving hard copies (optional printer)



Sophisticated data analysis capabilities with RPCR software

Cloud based paperless software

T-Top devices feature the exclusive and technologically advanced T-Connect, a cloud based IoT solution for your sterilization center.



Redefining Your Sterilization Center And Smart Clinics

T-Connect is a cloud based IoT solution designed for your sterilization center. The software provides clinics with data intelligence in real-time, accessible from anywhere, anytime, in a secure way.

Connect autoclaves, washers and more with Wi-Fi or the T-Connect gateway, and use the advanced T-Connect software to improve workflow efficiency and results.


Benefit from paperless documentation, device diagnostics, monitoring and troubleshooting. Knowing where and what the technical glitches are, allows you to have accelerated service, minimizing your device downtime.

Why Use T-Connect?

The ability to connect devices in your sterilization center and communicate on-site or remotely. Using the T-Connect gateway solution, it will connect to any device that is not Wi-Fi enabled.



What will T-Connect Do?



Connect To The Future Of Sterilization With T-Connect

- ◆ Connectivity by PC, tablet and smartphone
- ◆ Remotely access and download cycle history, get notifications and alerts
- ◆ Technicians and service centers can remotely analyze autoclaves, diagnose and resolve problems
- ◆ Suggestions on how to resolve the problem (video/instructions)
- ◆ Automatic notification to replenish consumables, remote diagnostics
- ◆ Data is safe
- ◆ Easy to use